



Rondish

SafeTime & SilverLining Pressure Sensor Pads/Mats


NEW!

Exciting, new-generation pressure pad sensors are now available - applying the latest technology for a safe, hygienic and very hard-wearing product. Working with our extensive range of electronic monitors, they are suitable for use in institutions or the homes, and carry a warranty of up to 1 year.

We offer a wide range of bed and chair pads, which sense when weight is removed or, alternatively, floor mats which sense when weight is put on. All pad sensors connect to our wired and wireless Bed/ChairWatcher alarm monitors, from standard carbon-based short-term (90 day) pads, to SilverLining long-life (1 year) pads and our latest, exciting innovation: SafeTime. SafeTime means caregivers will be warned when the 90-day or 1-year pads have exceeded their warranted lifetime, so there is no chance that a pad is forgotten about and deteriorates until it is in danger of failing. This innovation is customer-driven, based on feedback by caregivers in the field. It is effectively an electronic "use-by" alarm, ensuring healthcare and caregiving professionals can virtually guarantee the operational integrity and safety of their services, giving them confidence they will be alerted to come to the aid of vulnerable patients at all times.




STPBM-01 (30"x10")



BPL-01 (30"x20")



STPBS3-01 (15"x10")



FMAT-02 (24"x48")

SafeTime & SilverLining Pressure Sensor Pads/Mats

- SafeTime provides optimum reliability for caregivers and patients, automatically warning caregivers when pad life has expired. When timed out, SafeTime pads continue to work but emit a continuous, gentle beep with the "change pad" LED on.
- SafeTime pads have an internal clock which begins to time the life of the pad as soon as its plugged into a monitor. It therefore only counts when in use.
- SilverLining provides ultra-sensitive contact and a very hard-wearing surface. The result is virtually no false alarms over a long life.
- Our Bed and Chair pads use highly anti-allergenic soft PVC for patient comfort.
- Pads have vent holes in the neck away from the main area to prevent water ingress, plus strain-relief grip on the cable attachment.
- Instructions for care and hygiene are printed on the pads.
- Floor mats are pads that sense when weight is put on, for use next to beds or doors to prevent wandering and may also be used with wired and wireless Bed/ChairWatcher monitor alarms. Made of heavy-duty, non-slip material, designed to lie flat and last.
- Many wireless Bed/ChairWatcher monitors may be monitored remotely from the caregiver's desk using a Central Station.



Specifications

Model No.	Description	Size	Warranty period
BPL-01	SilverLining™ Pressure Sensor Bed Pad	30"x20"	1 year
BPM-01	SilverLining™ Pressure Sensor Bed Pad	30"x10"	1 year
BPS-01	SilverLining™ Pressure Sensor Wheelchair Pad	15"x10"	1 year
FMAT-01	SilverLining™ Light-weight Floor Mat.	30"x20"	1 year
FMAT-02	SilverLining™ Heavy-duty non-slip Floor Mat.	48"x24"	1 year
TP-01	Toilet Seat Pad	5.5"x1.5"	30 days
BP-01	Pressure Sensor Bed Pad	30"x20"	90 days
WP-01	Pressure Sensor wheelchair Pad	15"x10"	90 days
BP6-01	Pressure Sensor Bed Pad	30"x10"	6 months
STPBL-01	SafeTime & SilverLining™ Pressure Sensor Bed Pad	30"x20"	1 year
STPBM-01	SafeTime & SilverLining™ Pressure Sensor Bed Pad	30"x10"	1 year
STPBS-01	SafeTime & SilverLining™ Pressure Sensor Wheelchair Pad	15"x10"	1 year
STPBL3-01	SafeTime Pressure Sensor Bed Pad	30"x20"	90 days
STPBM3-01	SafeTime Pressure Sensor Bed Pad	30"x10"	90 days
STPBS3-01	SafeTime Pressure Wheelchair Pad	15"x10"	90 days

Note: The default lifetimes (90-day and 1-year as used in North America) are shown in the table but customer specification for batch orders of 1000 or more for each type of pad. Floor mats minimum order is 1500

Distributed By:



Hong Kong Office:
Rondish Company Limited
Tel: +852 2545 2057
Fax: +852 2541 7411
info@rondish.com
www.rondish.com



USA Office:
RondishUSA LLC
Tel: +1 858 729 0251
Fax: +1 815 346 5396
sales@rondishUSA.com
www.rondishUSA.com